
THUASNE

Stomex

fr	Ceinture abdominale pour stomie	6
en	Abdominal support for stoma	7
de	Bauchgurt für Stoma	8
nl	Buikgordel voor stoma	9
it	Cintura addominale per stomizzati	10
es	Faja abdominal para estoma	12
pt	Cinta abdominal para ostomia	13
da	Mavebælte til stomi	14
fi	Tukivyö vatsa-avanteelle	15
sv	Bukbälte för stomi	16
el	Ζώνη κοιλιακής χώρας για κολοστομία	18
cs	Stomický břišní pás	19
pl	Pas brzuszny do stomii	20
lv	Vēdera daļas josta stomas pacientiem	21
lt	Pilvo diržas esant ostomijai	22
et	Kõhupiirkonna tugivöö stoomi korral	23
sl	Oporni pas za trebušni predel in stomo	25
sk	Brušný podporný pás pre stomických pacientov	26
hu	Hasi támasztó öv sztómás betegeknek	27
bg	Колан за коремната област при стоматит	28
ro	Centură abdominală pentru stomă	29
ru	Корсет для брюшного отдела с отверстием для стомы	30
hr	Abdominalni pojas za stomu	32
zh	用于气孔的腹部支撑带	33
ar	حزام فغر البطن	35

	1	60 - 70 cm
	2	71 - 85 cm
	3	86 - 100 cm
	4	101 - 120 cm

fr	Tour de taille*	Hauteurs : 18 cm - 25 cm
en	Waist measurement*	Heights: 18 cm - 25 cm
de	Taillenumfang*	Höhen: 18 cm - 25 cm
nl	Tailleomvang*	Hoogtes: 18 cm - 25 cm
it	Circonferenza vita*	Altezze: 18 cm - 25 cm
es	Contorno de cintura*	Alturas: 18cm - 25cm
pt	Perímetro de cintura*	Alturas: 18cm - 25cm
da	Taljemål*	Højde: 18cm - 25cm
fi	Vyötärönympärys*	Korkeudet: 18 cm - 25 cm
sv	Midjemått*	Höjd: 18cm - 25cm
el	Περίμετρος μέσης*	Ύψη: 18 εκ. - 25 εκ.
cs	Obvod pasu*	Výšky: 18 cm - 25 cm
pl	Obwód w talii*	Wysokości: 18 cm - 25 cm
lv	Vidukļa apkārtmērs*	Augstumi: 18 cm - 25 cm
lt	Juosmens apimtis*	Aukštis: 18 cm - 25 cm
et	Vööümbermõõt*	Kõrgused: 18 cm - 25 cm
sl	Obseg pasu*	Višina: 18 cm - 25 cm
sk	Obvod pásu*	Výšky: 18 cm - 25 cm
hu	Derékbőség*	Magasság: 18 cm - 25 cm
bg	Обиколка на талията*	Височини: 18 cm - 25 cm
ro	Circumferința taliei*	Înălțimi: 18 cm - 25 cm
ru	Обхват талии*	Высота: 18 см - 25 см
hr	Opseg struka*	Visine: 18 cm - 25 cm
zh	腰围*	高度 : 18 cm-25 cm
ar	محيط الخصر*	الارتفاعات: 18 سم - 25 سم

Compression

Soutien

Compression

Support

Kompression

Stütze

Compressie

Ondersteuning

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Sostegno

Compresión

Sujeción

Compressão

Suporte

Kompression

Støtte

Kompressio

Tuki

Kompression

Stöd

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ABDOMINAL SUPPORT FOR STOMA

Description/Destination

The device is intended only for the treatment of the indications listed and for patients whose measurements correspond to the sizing table*.

Composition

Compact elastic fabric. ①

Abdominal panel in elastic, ladder-proof fabric to enable openings to be cut out. ②

Closure with self-fastening surface and thumb-loops. ③

Flexible and removable dorsal stays. ④

Flexible and removable lateral stays. ⑤

Textile components: cotton - elastane - elastodiene - polyamide - polyester - polyether - viscose.

Non-textile components: brushed tempered steel - polyamide

Properties/Mode of action

The device ensures strong support and adaptable compression to the abdomen thanks to the elastic compact fabric.

Indications

Stoma.

Post-operative abdominal conditions.

Contraindications

Do not apply the product in direct contact with broken skin.

Do not use in the event of known allergy to any of the components.

Do not use on pregnant women.

Do not use in the event of a hiatus hernia.

Do not use in the event of bone cancer with spinal metastases.

Do not use in the event of circulatory, pulmonary or cardiovascular problems in patients for whom an increase in arterial blood pressure is not recommended.

Precautions

Verify the product's integrity before every use.

Do not use the device if it is damaged.

Choose the appropriate size to fit the patient, referring to the size chart.

Strictly comply with your healthcare professional's prescription and recommendations for use.

It is recommended that a healthcare professional adjust the stays to the shape of the patient's back before the belt is worn for the first time.

It is recommended that a healthcare professional supervises the first application.

Do not wear the product in a medical imaging machine.

Do not wear the product while sleeping.

In the event of discomfort, significant hindrance, pain, abnormal sensations, remove the device and consult a healthcare professional.

For hygiene and performance reasons, do not re-use the product for another patient.

It is recommended to adequately tighten the device to achieve support without excessive compression.

Do not use the device in case of application of certain products on the skin (creams, ointments, oils, gels, patches...).

The use of a medical device by a child, should be done under the supervision of an adult or a healthcare professional.

Undesirable side-effects

This device can cause skin reactions (redness, itching, burns, blisters, etc.) or wounds of various degrees of severity.

In postoperative conditions, the restrictive effect of the device could alter respiratory function and decrease lung capacity and tidal volume.

Any serious incidents occurring related to the device should be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is resident.

Instructions for use/Application

Fit the belt so as to determine where to cut the opening.

Position the belt around the waist, using the thumb-loops. ④

Outline the opening on the abdominal panel of the belt at the level of the stoma. ⑥

Remove the belt and cut out with scissors the opening. ⑦

Fitting the device

Position the belt around the waist, using the thumb-loops.

The belt must be centred on the spine: the dorsal stays should be positioned on either side of the spine.

The bottom of the belt must stop at the top of the crease between the buttocks.

Hold the left end flat using your right hand.

Pull the thumb-loops simultaneously and evenly. Ⓞ

Fasten the belt, left side over right side, adjusting to the desired compression. Ⓞ

If the self-fastening strap is properly positioned it does not wrinkle.

The level of compression can be adjusted as required throughout the day.

Care

Close the self-fastening tabs before washing. Machine washable at 60 °C (delicate programme). Do not machine wash the device more than 10 times. Do not dry clean. Do not use detergents, fabric softeners or aggressive products (products containing chlorine). Do not tumble-dry. Do not iron. Squeeze out excess water. Dry away from any direct heat source (radiator, sun, etc.). Dry flat. If the belt gets dirty, it is recommended and, exceptionally, possible to wash it at 90 °C.

Storage

Store at room temperature, preferably in the original packaging.

Disposal

Dispose of in accordance with local regulations.

* Not applicable to made-to-measure medical devices, where measurements will be taken individually for each patient.

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Nusiimkite diržą ir žirkėmis iškirpkite angą. ©