THUASNE

PRO MASTER

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Dear patient

THUASNE has developed the innovative **PROMASTER** bandages as modern orthopaedic aids to help you lead an active and mobile life without discomfort. All therapeutic products in the **PROMASTER** range set standards in respect of shape, styling and functional performance. They are exceedingly comfortable and provide the reassurance you need, even when exercising.

The **Genu PROMASTER** bandage promotes the best possible healing of your injured or painful knee. Place your trust in this premium bandage featuring the very latest technology.

Genu PROMASTER

Stimulates the muscles. Relieves pressure points.

The **Genu PROMASTER** knee bandage features an adaptable Y-shaped silicone pad. It was designed to particularly stimulate proprioception – the unconscious perception of movements in the area of the thigh muscle insertions, which speeds up the healing process. The **Genu PROMASTER** gently corrects posture and stabilises the knee.

- Fast and lasting pain reduction
- Y pad relieves tensions and supports the patella
- FlexAir pad fits perfectly and transports moisture away from the skin

Comfortable, firm fit

- Gradual change in compression prevents pinching
- Breathable knit fabric with 3-D texture for comfort

Y pad

- Patella support
- Surrounding tissue is massaged, extending right up to the thighs
- The lower taper of the pad is cut out to relieve pressure on the edge of the shinbone

Spiral springs

- Enhanced stabilisation of the knee joint
- Gently support knee movement without obstructing it

Pull-on aids on both sides

- Make it easier to apply the bandage

Exercise regime

Knee

We have developed a short exercise regime for you to promote the healing process. The exercises are designed so that you can perform them anywhere, any time with normal everyday aids. They will take between five and ten minutes, but you can extend or shorten them by altering the number of repetitions. Regular exercise will restore the stability and functional reliability of your knee more quickly.

EXERCISE AIDS

Towel

Chair

Cushion

NOTES

- Before starting the exercises, make sure you are standing securely.
- Warm up before exercising: don't stretch cold muscles!
- Stretch gently and carefully just until you begin to feel pain; never go beyond the pain threshold.
- Perform the exercises slowly; focus on what you are doing.
- Continue to breathe calmly and evenly during the exercises.

CAUTION!

None of the exercises should cause pain or exacerbate existing pain. If you have any questions about the exercises themselves, the intensity or frequency, etc., please consult your doctor or physiotherapist.


EXERCISE 1

Warm up – Walking on the spot

Lift your heels one after the other from the floor. Make sure your knees move forward during the lift, not to the side. Try to keep the movements symmetrical on both sides. To increase the intensity, try standing on a cushion to exercise.

Duration:

Walk for 60 seconds – pause for 30 seconds
– walk for 60 seconds


EXERCISE 2

Stretching the knee extensor


Bend your supporting leg slightly. Pull your heel towards your bottom using either your hand or a towel. Keep your back straight, avoid a hollow back.

Duration:

Hold for 2 x 15 seconds on each side.

EXERCISE 3


Stretching the rear knee joint

Step back with the affected leg, keeping your foot on the ground. Extend your knee joint until you feel the muscles stretching slightly in the front leg.

Duration:

Hold for 2 x 15 seconds on each side.

EXERCISE 4

Strengthening - Squats

Stand up straight with your feet shoulder-width apart. Hold onto the back of a chair for stability. The chair is only meant to help you keep your balance. Slowly squat, but only as far as you can without feeling pain. Look straight ahead and slightly up. Make sure you don't hollow your back while doing the exercises. Don't squat too low to start with.

Duration:

3 x 10 repetitions

EXERCISE 5


Balance and strengthening the postural muscles

Stand on the afflicted leg. At the start of the exercise, hold onto the back of a chair for safety. Then lift the healthy leg and bend the knee. Slowly swing the leg back and forth. Keep standing upright and keep the supporting leg straight. Repeat the exercise with the other leg.

Try to do the exercise without holding on to the chair, but only if you feel confident enough.

Duration:

30 repetitions on each side

We hope your knee heals quickly and well.


Your **THUASNE** Team

MEIN PERSÖNLICHER TRAININGSPLAN MY PERSONAL TRAINING SCHEDULE

Übungen/Exercises:

.....X Wiederholungen/Sets des
Komplettprogramms
Repetitions/sets of the
complete exercise regime

.....X pro Woche/per week:

Meine Trainingstage/my training days

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THUASNE DEUTSCHLAND GmbH