THUASNE®

Sleeq Flex

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Nylon, Polyester, Polypropylen, synthetischer Kautschuk (latexfrei)

Sleeq Flex (54 660)

Spinal orthosis for active relief and correction of the lumbar and the thoracic spine in the sagittal plane

Purpose

The Sleeq Flex spinal orthosis is an orthosis for active relief and correction of the lumbar and the thoracic spine in the sagittal plane. The orthosis may only be used for lower spine treatment.

Properties

- Pulley system for individual power transmission
- Pliable material
- Reinforcing elements in the abdominal area

Indications

- Stable, osteoporotic vertebral fractures
- Hunchback with poor shoulder posture
- Spondylitis
- Morbus Bechterew

Contraindications

- Unstable vertebral fractures
- Do not apply the product near open wounds
- Haematomas and severe swelling

Side effects

When properly applied, no side effects have been known to date.

Putting on the orthosis

Size determination and adjustment

The Sleeq Flex spinal orthosis has a universal size and can be individually adjusted to the circumferences.

(Fig. 1) Four size markings (1, 2, 3, 4) show markings for the waist size and the correct alignment. They are merely used for orientation since the size can be varied at will and individually adjusted.

(Fig. 2) Loosen the tension straps and place them on the spinal pad, with the hook-and-loop fastening facing up. Now push the right wing to the inside or the outside and firmly press the crescent-shaped hook-and-loop fastening on the suitable size marking (or between the markings). Repeat this setting with the left wing as well. Make sure that the hook-and-loop fastening fits firmly between the waist strap and the crescent-shaped hook-and-loop part. If the desired setting does not fit, repeat it until the circumference is correct.

(Fig. 3) If necessary, the waist strap may be fixed slightly diagonally upwards or downwards to thus obtain a better waist adjustment.

In order to relieve the pulley system, leave the tension straps still on the spinal pad and simultaneously pull on the left and right waist strap. Once you feel any resistance, the system will be adjusted and both waist straps should now have the same length. You can now firmly fasten, by hook and loop, the tension straps on the left and right each (without pull) on the waist strap.

Putting on the orthosis

(Fig. 4) Loosen the shoulder straps and push the arms through the straps (for larger patients, the „ThoraCare“ carrier may be pushed upwards before putting it on). Now centre the spinal pad at the desired height in the lumbar region and then close the ventral wings (without pull). To this end, place the left strap on the abdomen and use the

en

hook-and-loop system to fasten the right strap on the left.

(Fig. 5a) For the individual pressure setting, please loosen the tension straps on the outside. Then put your thumbs into the opening and pull the straps forward to the outside until you feel a comfortable and firming pressure. In order to fix the two strap ends, simply use the hook-and-loop system to fasten the flaps again on the waist strap. In order to firm up the shoulder straps, simply pull the traction ropes downwards (backpack system) until they rest comfortably against the shoulders.

(Fig. 5b) In order to adjust the rear lengths, push the „ThoraCare“ carrier upward or downward in the pocket so that the upper centre of the carrier is on a level with T3.

Precautions

Follow the advice of the professional who prescribed or supplied the product. Consult your professional if discomfort occurs. Store at room temperature, preferably in its original packaging. For hygiene and performance reasons, do not re-use the product for another patient. Any serious incident that occurs in connection with the product must be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient are/is domiciled.

Care

Do not use any cleaning agents, conditioners or aggressive products (containing chlorine). Squeeze out the water thoroughly. Dry away from any heat sources

Material composition

Nylon, polyester, polypropylene, synthetic rubber (latex-free)

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Thuasne Cervitex Ltd. · The Bell House, 20 Hataas St. · Kfar Saba, 44425
Thuasne Scandinavia AB · Box 95 · 131 07 Nacka
Thuasne UK Ltd
Thuasne Deutschland GmbH