**THUASNE®**

# PRO MASTER

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Ihr **Thuasne**® Team

Dear patient

Thusne® has developed the innovative **PROMASTER** bandages as modern orthopaedic aids to help you lead an active and mobile life without discomfort. All therapeutic products in the **PROMASTER** range set standards in respect of shape, styling and functional performance. They are exceedingly comfortable and provide the reassurance you need, even when exercising.

The **Malleo PROMASTER** bandage promotes the best possible healing of your injured or painful ankle. Place your trust in this premium bandage featuring the very latest technology.

# Malleo PROMASTER

## Firm hold. All-round support.

The anatomically shaped **Malleo PROMASTER** supports the ankle on both sides. Even compression, supported by point-elastic pads, reduces oedema and pain. Soft, highly flexible knit fabric for enhanced comfort.

- Air-permeable, point-elastic silicone pads reduce oedema
- Breathable fabric
- Decreased compression at the edges of the bandage prevent pinching

**Softly tapered change  
in compression**

- Prevents pinching

**Pull-on aids  
on both sides**

- Make it easier to put the bandage on

**Anatomically shaped  
pads**

- Effectively support the ankle

**Pronation straps (only on Malleo PROMASTER Aktiv)**

- Enhanced forefoot stability and protection against the foot buckling inwards
- Thus preventing resulting additional stress on the ankle ligaments and tendons

# Exercise regime

## Ankle

We have developed a short exercise regime for you to promote the healing process. The exercises are designed so that you can perform them anywhere, any time with normal everyday aids. They will take between five and ten minutes, but you can extend or shorten them by altering the number of repetitions. Regular exercise will restore the stability and functional reliability of your ankle more quickly.

### EXERCISE AIDS

Chair

Cushion


### NOTES

- Before starting the exercises, make sure you are standing or sitting securely.
- Warm up before exercising; don't stretch cold muscles!
- Stretch gently and carefully just until you begin to feel pain; never go beyond the pain threshold.
- Perform the exercises slowly; focus on what you are doing.
- Continue to breathe calmly and evenly during the exercises.

### CAUTION!

None of the exercises should cause pain or exacerbate existing pain. If you have any questions about the exercises themselves, the intensity or frequency, etc., please consult your doctor or physiotherapist.

## EXERCISE 1


### Warm up – Walking on the spot

Lift your heels one after the other from the floor. Make sure your knees move forward during the lift, not to the side. Try to keep the movements symmetrical on both sides.

**Duration:**

Walk for 60 seconds – pause for 30 seconds  
– walk for 60 seconds

## EXERCISE 2

### Stretching the calf muscles

Step back with the affected leg, keeping the sole of your foot on the ground. Extend the knee of your back leg.

**Duration:**

Hold for 2 x 20 seconds on each side

### EXERCISE 3


#### Strengthening the foot lifting muscles

Sit on a chair with your knees bent at right angles and your feet flat on the floor. Bend both ankles up as far as they will go, with the toes pointing upwards and the heels on the ground.

Hold for approx. 5 seconds, then release for 5 seconds.

**Duration:**

5 repetitions, more if necessary

### EXERCISE 4

#### Balance and strengthening the postural muscles


Stand on the afflicted foot.

At the start of the exercise, hold onto the back of a chair for safety. Then lift the healthy leg and bend the knee. Slowly swing the leg back and forth. Keep standing upright and keep the supporting leg straight.

Try to do the exercise without holding on to the chair, but only if you feel confident enough.

**Duration:**

30 repetitions on each side

## EXERCISE 5



### Coordination – Stepping sequence

Lift your heels one after the other from the floor, rolling along the outer sides of your feet.

Try to keep the movements symmetrical on both sides. To increase the intensity, try standing on a cushion to exercise.

**Duration:**

30 repetitions on each side

We hope your ankle heals quickly and well.


Your **Thuasne**® Team

# MEIN PERSÖNLICHER TRAININGSPLAN MY PERSONAL TRAINING SCHEDULE

Übungen/Exercises:

.....X Wiederholungen/Sets des  
Komplettprogramms  
Repetitions/sets of the  
complete exercise regime

.....X pro Woche/per week:

Meine Trainingstage/my training days

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**Thuasne Deutschland GmbH**