THUASNE®

PRO
MASTER

Genu PFS PROMASTER
Genu FS

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
Ihr **Thuasne**® Team

Dear patient

Thusne® has developed these innovative bandages as modern orthopaedic aids to help you lead an active and mobile life.

These therapeutic products set standards in respect of shape, styling and functional performance. They are exceedingly comfortable and provide the reassurance you need, even when exercising.

The Genu PFS PROMASTER and Genu FS orthosis support the stability and functional reliability of your knee. Place your trust in these premium orthoses featuring the very latest technology.

A close-up photograph of a dark, textured surface, possibly a shoe or bag, with a prominent blue diagonal stripe in the upper right corner. The texture is a fine, repeating pattern of small, raised, rounded shapes.

Genu PFS PROMASTER
and Genu FS
Reliable support.
Specific correction.

Exercise regime

Knee

The following exercises will improve the stability and functional reliability of your knee. The exercises are designed so that you can perform them anywhere, any time with normal everyday aids. They will take between five and ten minutes, but you can extend or shorten them by altering the number of repetitions. Regular exercise will restore the stability and functional reliability of your knee more quickly.

EXERCISE AIDS

Towel

Chair

2 cushions


NOTES

- Before starting the exercises, make sure you are standing securely.
- Warm up before exercising: don't stretch cold muscles!
- Stretch gently and carefully just until you begin to feel pain; never go beyond the pain threshold.
- Perform the exercises slowly; focus on what you are doing.
- Continue to breathe calmly and evenly during the exercises.

CAUTION!

None of the exercises should cause pain or exacerbate existing pain. If you have any questions about the exercises themselves, the intensity or frequency, etc., please consult your doctor or physiotherapist.

EXERCISE 1


Warm up – Walking on the spot

Lift your heels one after the other from the floor. Make sure your knees move forward during the lift, not to the side. Try to keep the movements symmetrical on both sides. To increase the intensity, try standing on a cushion to exercise.

Duration:

Walk for 60 seconds – pause for 30 seconds – walk for 60 seconds

EXERCISE 2


Stretching the knee extensor

Bend your supporting leg slightly. Pull your heel towards your bottom using either your hand or a towel. Keep your back straight, avoid a hollow back.

Duration:

Hold for 2 x 15 seconds on each side.

EXERCISE 3

Stretching the rear knee joint

Step back with the affected leg, keeping your foot on the ground. Extend your knee joint until you feel the muscles stretching slightly in the front leg.

Duration:

Hold for 2 x 15 seconds on each side.

EXERCISE 4

Squats with cushion

Stand up straight. Squeeze 2 small cushions between your thighs and your knees. Hold onto the back of a chair for stability. The chair is only meant to help you keep your balance. Slowly squat, but only as far as you can without feeling pain. Keep looking straight ahead and slightly up

to make sure your head remains straight. Make sure you don't hollow your back while doing the exercises. Don't squat too low to start with.

Duration: 3 x 10 repetitions

EXERCISE 5


Strengthening the inner thigh muscles

Squeeze a towel or cushion between your knees and press your knees and lower legs together firmly.

Hold for 5-10 seconds before releasing briefly.

Duration: 5 repetitions

EXERCISE 6

Balance and strengthening the postural muscles

Stand on the afflicted leg. At the start of the exercise, hold onto the back of a chair for safety. Then lift the healthy leg and bend the knee. Slowly swing the leg back and forth. Keep standing upright and keep the supporting leg straight. Repeat the exercise with the other leg.

Try to do the exercise without holding on to the chair, but only if you feel confident enough.

Duration:

30 repetitions on each side

EXERCISE 7

Strengthening the functional performance of the inner thigh muscles

Step to the side and back with your afflicted leg. Pull your toes up and point the tips of your toes outwards. Then slowly guide the leg to the front and centre (diagonal movement) and back to the starting position.

Duration:

5 repetitions on each side

EXERCISE 8

Strengthening the functional performance of the inner thigh muscles

Bend your healthy knee and place your foot on the floor. Stretch the afflicted leg out to the side with the toes pointing outwards. Swing the stretched leg over the healthy leg, as if you were at a small hill.

Duration: 5 repetitions on each side

We hope you enjoy the exercises and that they help your recovery.


Your **Thuasne**[®] Team

MEIN PERSÖNLICHER TRAININGSPLAN MY PERSONAL TRAINING SCHEDULE

Übungen/Exercises:

.....X Wiederholungen/Sets des
Komplettprogramms
Repetitions/sets of the
complete exercise regime

.....X pro Woche/per week:

Meine Trainingstage/my training days

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Thuasne Deutschland GmbH